





STPG2 Antibody

Product Code	CSB-PA847628LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N412
Immunogen	Recombinant Human Sperm-tail PG-rich repeat-containing protein 2 protein (201-459AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sperm-tail PG-rich repeat-containing protein 2, STPG2, C4orf37
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	STPG2

Image



Western blot
 All lanes: STPG2 antibody at 0.8µg/ml + HL60 whole cell lysate
 Secondary
 Goat polyclonal to rabbit at 1/10000 dilution
 Predicted band size: 51 kDa
 Observed band size: 51 kDa



Immunohistochemistry of paraffin-embedded human spleen tissue using CSB-PA847628LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.