



# VN1R2 Antibody, HRP conjugated

|                     |                                                                                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA847711LB01HU                                                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                   |
| Uniprot No.         | Q8NFZ6                                                                                                          |
| Immunogen           | Recombinant Human Vomeronasal type-1 receptor 2 protein (115-170AA)                                             |
| Raised In           | Rabbit                                                                                                          |
| Species Reactivity  | Human                                                                                                           |
| Tested Applications | ELISA                                                                                                           |
| Form                | Liquid                                                                                                          |
| Conjugate           | HRP                                                                                                             |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                |
| Purification Method | >95%, Protein G purified                                                                                        |
| Isotype             | IgG                                                                                                             |
| Clonality           | Polyclonal                                                                                                      |
| Alias               | Vomeronasal type-1 receptor 2 (G-protein coupled receptor GPCR25) (hGPCR25) (V1r-like receptor 2), VN1R2, V1RL2 |
| Immunogen Species   | Homo sapiens (Human)                                                                                            |
| Research Area       | Neuroscience                                                                                                    |
| Target Names        | VN1R2                                                                                                           |