



VN1R2 Antibody, Biotin conjugated

Product Code	CSB-PA847711LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFZ6
Immunogen	Recombinant Human Vomeronasal type-1 receptor 2 protein (115-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vomeronasal type-1 receptor 2 (G-protein coupled receptor GPCR25) (hGPCR25) (V1r-like receptor 2), VN1R2, V1RL2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	VN1R2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.