



WDR66 Antibody, HRP conjugated

Product Code	CSB-PA848398LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBY9
Immunogen	Recombinant Human WD repeat-containing protein 66 protein (1-1154AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cilia- and flagella-associated protein 251 (CFAP251) (WD repeat-containing protein 66), WDR66, CFAP251
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	WDR66
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.