



LGALS14 Antibody, Biotin conjugated

Product Code	CSB-PA848400LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCE9
Immunogen	Recombinant Human Placental protein 13-like protein (31-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Placental protein 13-like (Charcot-Leyden crystal protein 2) (CLC2) (Galectin-14) (Gal-14), LGALS14, PPL13
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LGALS14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.