



FBXO25 Antibody, FITC conjugated

Product Code	CSB-PA848401LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCJ0
Immunogen	Recombinant Human F-box only protein 25 protein (1-291AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box only protein 25, FBXO25, FBX25
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FBXO25
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.