



GPX8 Antibody, Biotin conjugated

Product Code	CSB-PA848417LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TED1
Immunogen	Recombinant Human Probable glutathione peroxidase 8 protein (41-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable glutathione peroxidase 8 (GPx-8) (GSHPx-8) (EC 1.11.1.9), GPX8
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	GPX8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.