



GTSF1 Antibody, Biotin conjugated

Product Code	CSB-PA848827LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WW33
Immunogen	Recombinant Human Gametocyte-specific factor 1 protein (1-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gametocyte-specific factor 1 (Protein FAM112B), GTSF1, FAM112B
Immunogen Species	Homo sapiens (Human)
Research Area	Developmental Biology
Target Names	GTSF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.