



GCDH Antibody, Biotin conjugated

Product Code	CSB-PA849798LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92947
Immunogen	Recombinant Human Glutaryl-CoA dehydrogenase, mitochondrial protein (45-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutaryl-CoA dehydrogenase, mitochondrial (GCD) (EC 1.3.8.6), GCDH
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	GCDH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.