





PPP1R14A Antibody

Product Code	CSB-PA850252LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A00
Immunogen	Recombinant Human Protein phosphatase 1 regulatory subunit 14A protein (1-147AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase 1 regulatory subunit 14A (17 kDa PKC-potentiated inhibitory protein of PP1) (Protein kinase C-potentiated inhibitor protein of 17 kDa) (CPI-17), PPP1R14A, CPI17 PPP1INL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PPP1R14A

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA850252LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA850252LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.