



ZNF274 Antibody, FITC conjugated

Product Code	CSB-PA850299LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GC6
Immunogen	Recombinant Human Neurotrophin receptor-interacting factor homolog protein (1-320AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurotrophin receptor-interacting factor homolog (Zinc finger protein 274) (Zinc finger protein HFB101) (Zinc finger protein with KRAB and SCAN domains 19) (Zinc finger protein zfp2) (Zf2), ZNF274, ZKSCAN19
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF274
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.