



MRGPRX3 Antibody, Biotin conjugated

Product Code	CSB-PA850382LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB0
Immunogen	Recombinant Human Mas-related G-protein coupled receptor member X3 protein (276-322AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mas-related G-protein coupled receptor member X3 (Sensory neuron-specific G-protein coupled receptor 1/2), MRGPRX3, MRGX3 SNSR1 SNSR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MRGPRX3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.