



RPRD1A Antibody, HRP conjugated

Product Code	CSB-PA850404LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96P16
Immunogen	Recombinant Human Regulation of nuclear pre-mRNA domain-containing protein 1A protein (2-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulation of nuclear pre-mRNA domain-containing protein 1A (Cyclin-dependent kinase inhibitor 2B-related protein) (p15INK4B-related protein), RPRD1A, P15RS
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RPRD1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.