





# FERD3L Antibody

<b>Product Code</b>	CSB-PA850424LA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q96RJ6
<b>Immunogen</b>	Recombinant Human Fer3-like protein (1-167AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Fer3-like protein (Basic helix-loop-helix protein N-twist) (Class A basic helix-loop-helix protein 31) (bHLHa31) (Nephew of atonal 3) (Neuronal twist), FERD3L, BHLHA31 NATO3 NTWIST
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	FERD3L

## Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA850424LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA850424LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)