





FERD3L Antibody

Product Code	CSB-PA850424LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RJ6
Immunogen	Recombinant Human Fer3-like protein (1-167AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fer3-like protein (Basic helix-loop-helix protein N-twist) (Class A basic helix-loop-helix protein 31) (bHLHa31) (Nephew of atonal 3) (Neuronal twist), FERD3L, BHLHA31 NATO3 NTWIST
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FERD3L

Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA850424LA01HU at dilution of 1:100



Immunofluorescent analysis of HeLa cells using CSB-PA850424LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.