




GALNT16 Antibody

Product Code	CSB-PA850796LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N428
Immunogen	Recombinant Human Polypeptide N-acetylgalactosaminyltransferase 16 protein (301-558AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GALNT16

Image



Western blot
All lanes: GALNT16 antibody at 1:500+ Mouse kidney tissue
Secondary
Goat polyclonal to rabbit at 1/10000 dilution
Predicted band size: 64, 61 kDa
Observed band size: 64, 49 kDa