



PRELID2 Antibody

Product Code	CSB-PA850836LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N945
Immunogen	Recombinant Human PRELI domain-containing protein 2 protein (1-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PRELI domain-containing protein 2, PRELID2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PRELID2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.