



EIF1AD Antibody, HRP conjugated

Product Code	CSB-PA850842LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9N8
Immunogen	Recombinant Human Probable RNA-binding protein EIF1AD protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable RNA-binding protein EIF1AD (Eukaryotic translation initiation factor 1A domain-containing protein) (Haponin), EIF1AD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EIF1AD