



EIF1AD Antibody, Biotin conjugated

Product Code	CSB-PA850842LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9N8
Immunogen	Recombinant Human Probable RNA-binding protein EIF1AD protein (1-165AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable RNA-binding protein EIF1AD (Eukaryotic translation initiation factor 1A domain-containing protein) (Haponin), EIF1AD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EIF1AD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.