



TDRD9 Antibody, Biotin conjugated

Product Code	CSB-PA850866LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDG6
Immunogen	Recombinant Human Putative ATP-dependent RNA helicase TDRD9 protein (286-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase TDRD9 (EC 3.6.4.13) (Tudor domain-containing protein 9), TDRD9, C14orf75
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TDRD9