



TDRD9 Antibody, Biotin conjugated

Product Code	CSB-PA850866LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDG6
Immunogen	Recombinant Human Putative ATP-dependent RNA helicase TDRD9 protein (286-585AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-dependent RNA helicase TDRD9 (EC 3.6.4.13) (Tudor domain-containing protein 9), TDRD9, C14orf75
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TDRD9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.