



NFAM1 Antibody, HRP conjugated

Product Code	CSB-PA850874LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NET5
Immunogen	Recombinant Human NFAT activation molecule 1 protein (186-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NFAT activation molecule 1 (Calcineurin/NFAT-activating ITAM-containing protein) (NFAT-activating protein with ITAM motif 1), NFAM1, CNAIP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NFAM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.