



NFAM1 Antibody, HRP conjugated

Product Code	CSB-PA850874LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NET5
Immunogen	Recombinant Human NFAT activation molecule 1 protein (186-270AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NFAT activation molecule 1 (Calcineurin/NFAT-activating ITAM-containing protein) (NFAT-activating protein with ITAM motif 1), NFAM1, CNAIP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NFAM1