



sodB Antibody

Product Code	CSB-PA852873LA01NGS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92450
Immunogen	Recombinant Neosartorya fumigata Superoxide dismutase [Mn], mitochondrial protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Neosartorya fumigata
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Superoxide dismutase [Mn], mitochondrial (EC 1.15.1.1) (allergen Asp f 6), sodB
Immunogen Species	Neosartorya fumigata (strain ATCC MYA-4609 / Af293 / CBS 101355 / FGSC A1100) (Aspergillus fumigatus)
Research Area	Others
Target Names	sodB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.