



RPE Antibody, Biotin conjugated

Product Code	CSB-PA853385LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AT9
Immunogen	Recombinant Human Ribulose-phosphate 3-epimerase protein (2-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribulose-phosphate 3-epimerase (EC 5.1.3.1) (Ribulose-5-phosphate-3-epimerase), RPE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPE
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.