



SARAF Antibody, Biotin conjugated

Product Code	CSB-PA853392HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BY9
Immunogen	Recombinant Human Store-operated calcium entry-associated regulatory factor protein (195-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Store-operated calcium entry-associated regulatory factor (SARAF) (SOCE-associated regulatory factor) (HBV X-transactivated gene 3 protein) (HBV XAg-transactivated protein 3) (Protein FOAP-7) (Transmembrane protein 66), SARAF, TMEM66 XTP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SARAF