



SARAF Antibody, Biotin conjugated

Product Code	CSB-PA853392HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96BY9
Immunogen	Recombinant Human Store-operated calcium entry-associated regulatory factor protein (195-339AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Store-operated calcium entry-associated regulatory factor (SARAF) (SOCE-associated regulatory factor) (HBV X-transactivated gene 3 protein) (HBV XAg-transactivated protein 3) (Protein FOAP-7) (Transmembrane protein 66), SARAF, TMEM66 XTP3
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SARAF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.