



S100A16 Antibody, Biotin conjugated

Product Code	CSB-PA853418LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96FQ6
Immunogen	Recombinant Human Protein S100-A16 protein (1-103AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A16 (Aging-associated gene 13 protein) (Protein S100-F) (S100 calcium-binding protein A16), S100A16, S100F
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	S100A16