



MRAP2 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA853421LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96G30 |
| Immunogen | Recombinant Human Melanocortin-2 receptor accessory protein 2 protein (66-205AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Melanocortin-2 receptor accessory protein 2 (MC2R accessory protein 2), MRAP2, C6orf117 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | MRAP2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |