





ACSBG1 Antibody

Product Code	CSB-PA853424LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GR2
Immunogen	Recombinant Human Long-chain-fatty-acid--CoA ligase ACSBG1 protein (351-600AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Long-chain-fatty-acid--CoA ligase ACSBG1 (EC 6.2.1.3) (Acyl-CoA synthetase bubblegum family member 1) (hBG1) (hsBG) (hsBGM) (Lipidosin), ACSBG1, BGM KIAA0631 LPD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACSBG1

Image



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA853424LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA853424LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.