



SEC16B Antibody, FITC conjugated

Product Code	CSB-PA853439LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JE7
Immunogen	Recombinant Human Protein transport protein Sec16B protein (5-304AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein transport protein Sec16B (Leucine zipper transcription regulator 2) (Regucalcin gene promoter region-related protein p117) (RGPR-p117) (SEC16 homolog B), SEC16B, KIAA1928 LZTR2 RGPR SEC16S
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SEC16B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.