




FAM71F1 Antibody

Product Code	CSB-PA853446LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KD3
Immunogen	Recombinant Human Protein FAM71F1 protein (121-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM71F1 (Protein FAM137A) (Testis development protein NYD-SP18), FAM71F1, FAM137A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM71F1

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA853446LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.