



FAM71F1 Antibody, FITC conjugated

Product Code	CSB-PA853446LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KD3
Immunogen	Recombinant Human Protein FAM71F1 protein (121-343AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM71F1 (Protein FAM137A) (Testis development protein NYD-SP18), FAM71F1, FAM137A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM71F1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.