



SLC26A8 Antibody, HRP conjugated

Product Code	CSB-PA853496LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RN1
Immunogen	Recombinant Human Testis anion transporter 1 protein (414-595AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis anion transporter 1 (Anion exchange transporter) (Solute carrier family 26 member 8), SLC26A8, TAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC26A8