



SLC26A8 Antibody, Biotin conjugated

Product Code	CSB-PA853496LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RN1
Immunogen	Recombinant Human Testis anion transporter 1 protein (414-595AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Testis anion transporter 1 (Anion exchange transporter) (Solute carrier family 26 member 8), SLC26A8, TAT1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC26A8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.