



FAM89B Antibody, Biotin conjugated

Product Code	CSB-PA854045LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N5H3
Immunogen	Recombinant Human Leucine repeat adapter protein 25 protein (1-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine repeat adapter protein 25, FAM89B, Lrap25
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FAM89B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.