



NANP Antibody, HRP conjugated

Product Code	CSB-PA855041LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TBE9
Immunogen	Recombinant Human N-acylneuraminate-9-phosphatase protein (1-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-acylneuraminate-9-phosphatase (EC 3.1.3.29) (Haloacid dehalogenase-like hydrolase domain-containing protein 4) (Neu5Ac-9-Pase), NANP, C20orf147 HDHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NANP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.