



NANP Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA855041LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8TBE9 |
| Immunogen | Recombinant Human N-acylneuraminate-9-phosphatase protein (1-248AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | N-acylneuraminate-9-phosphatase (EC 3.1.3.29) (Haloacid dehalogenase-like hydrolase domain-containing protein 4) (Neu5Ac-9-Pase), NANP, C20orf147 HDHD4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | NANP |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |