



SPPL2A Antibody, Biotin conjugated

Product Code	CSB-PA855050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCT8
Immunogen	Recombinant Human Signal peptide peptidase-like 2A protein (26-172AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptide peptidase-like 2A (SPP-like 2A) (SPPL2a) (EC 3.4.23.-) (Intramembrane protease 3) (IMP-3) (Presenilin-like protein 2), SPPL2A, IMP3 PSL2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SPPL2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.