



DUSP19 Antibody, Biotin conjugated

Product Code	CSB-PA855469LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WTR2
Immunogen	Recombinant Human Dual specificity protein phosphatase 19 protein (1-217AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 19 (EC 3.1.3.16) (EC 3.1.3.48) (Dual specificity phosphatase TS-DSP1) (Low molecular weight dual specificity phosphatase 3) (LMW-DSP3) (Protein phosphatase SKRP1) (Stress-activated protein kinase pathway-regulating phosphatase 1) (SAPK pathway-regulating phosphatase 1), DUSP19, DUSP17 LMWDSP3 SKRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.