




SLC7A3 Antibody

Product Code	CSB-PA855505LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WY07
Immunogen	Recombinant Human Cationic amino acid transporter 3 protein (428-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cationic amino acid transporter 3 (CAT-3) (CAT3) (Cationic amino acid transporter y+) (Solute carrier family 7 member 3), SLC7A3, ATRC3 CAT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC7A3

Image



Immunofluorescent analysis of MCF-7 cells using CSB-PA855505LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.