



SLC7A3 Antibody, HRP conjugated

Product Code	CSB-PA855505LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WY07
Immunogen	Recombinant Human Cationic amino acid transporter 3 protein (428-475AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cationic amino acid transporter 3 (CAT-3) (CAT3) (Cationic amino acid transporter y+) (Solute carrier family 7 member 3), SLC7A3, ATRC3 CAT3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SLC7A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.