



# SLC7A3 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA855505LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8WY07
<b>Immunogen</b>	Recombinant Human Cationic amino acid transporter 3 protein (428-475AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Cationic amino acid transporter 3 (CAT-3) (CAT3) (Cationic amino acid transporter y+) (Solute carrier family 7 member 3), SLC7A3, ATRC3 CAT3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	SLC7A3