



VP35 Antibody

Product Code	CSB-PA856262LA01RCT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91DE0
Immunogen	Recombinant Reston ebolavirus Polymerase cofactor VP35 protein (1-329AA)
Raised In	Rabbit
Species Reactivity	Reston ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polymerase cofactor VP35, VP35
Immunogen Species	Reston ebolavirus (strain Philippines-96) (REBOV) (Reston Ebola virus)
Research Area	Others
Target Names	VP35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.