



KRT35 Antibody, HRP conjugated

Product Code	CSB-PA856422LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92764
Immunogen	Recombinant Human Keratin, type I cuticular Ha5 protein (1-455AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Keratin, type I cuticular Ha5 (Hair keratin, type I Ha5) (Keratin-35) (K35), KRT35, HHA5 HKA5 KRTHA5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KRT35
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.