




MS4A6E Antibody

Product Code	CSB-PA856919LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DS6
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 6E protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 6E, MS4A6E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A6E

Image



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA856919LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.