



MS4A6E Antibody, FITC conjugated

Product Code	CSB-PA856919LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DS6
Immunogen	Recombinant Human Membrane-spanning 4-domains subfamily A member 6E protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-spanning 4-domains subfamily A member 6E, MS4A6E
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MS4A6E