



OMA1 Antibody, Biotin conjugated

Product Code	CSB-PA856923LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96E52
Immunogen	Recombinant Human Metalloendopeptidase OMA1, mitochondrial protein (22-194AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloendopeptidase OMA1, mitochondrial (EC 3.4.24.-) (Metalloprotease-related protein 1) (MPRP-1) (Overlapping with the m-AAA protease 1 homolog), OMA1, MPRP1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	OMA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.