



PRAC1 Antibody

Product Code	CSB-PA856965LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KF2
Immunogen	Recombinant Human Small nuclear protein PRAC1 protein (1-57AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small nuclear protein PRAC1 (Prostate cancer susceptibility candidate protein 1) (Prostate, rectum and colon expressed gene protein), PRAC1, C17orf92 PRAC
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PRAC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.