



LRRC3B Antibody, Biotin conjugated

Product Code	CSB-PA857000LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96PB8
Immunogen	Recombinant Human Leucine-rich repeat-containing protein 3B protein (34-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing protein 3B (Leucine-rich repeat protein LRP15), LRRC3B, LRP15
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LRRC3B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.